

PO21: Precision Oncology Course

Panel Analysis - GastroIntestinal Cancer Patient

Study Case - Gastrointestinal Cancer Patient


Tumor type: Patient with lung metastasis (gastrointestinal stromal tumor)

Sequencing: Illumina HiSeq2500

Type of data: Sequencing panel (paired).
Ion Ampliseq Comprehensive Cancer Panel
(409 genes)

Reference Genome: hg19

Samples: One tumor sample without
normal tissue



Gastrointestinal Stromal Tumors (GIST), from
Mangalore Institute of Oncology.

Study Case - Panel 2

Tumor type: Patient with lung metastasis (gastrointestinal stromal tumor)

Sequencing: Illumina HiSeq2500

Type of data: Sequencing panel (paired).
Ion Ampliseq Comprehensive Cancer Panel
(409 genes)

Reference Genome: hg19

Samples: One tumor sample without
normal tissue



hg19



fastq_panel2_tumor



tumor_1.fq



tumor_2.fq



bedfiles

Move to working directory

1. Create a new folder in de Desktop and move to it

```
$ mkdir Desktop/GIST
```

```
$ cd Desktop/GIST
```

2. Clone varca

```
$ git clone https://gitlab.com/bu\_cnio/varca
```

3. Move inside varca folder

```
$ cd varca
```

Edit configuration files

contigs.tsv

We are leaving this file empty so the analysis is performed using all contigs of the fasta index.

```
$ touch contigs.tsv
```

Edit configuration files

samples.tsv

As we only have tumor sample, MuTect2 will be executed in **tumor-only mode**. In order to do so, the **sample name** must be **the same in the sample and control columns**.

sample	control
tumor	tumor

Rename samples file: **\$ mv samples-example.tsv samples.tsv**

Edit configuration files

units.tsv

We only have **one tumor sample** with **two FASTQ files**: read 1 and read2.

sample	unit	platform	fq1	fq2
tumor	1	ILLUMINA	/Users/participant/Desktop/ panel2/tumor_1.fq	/Users/participant/Desktop/ panel2/tumor_2.fq

Rename units file: **\$ mv units-example.tsv units.tsv**

Edit configuration files

config.yaml

```
$ cp config-example.yaml config.yaml
```

```
$ vi config.yaml
```

```
samples: samples.tsv
```

```
units: units.tsv
```

```
contigs: contigs.tsv
```

```
outdir: "out"
```

```
logdir: "log"
```

```
ref:
```

```
# Genome database of snpeff to be used in the annotation with this resource. Available databases can be checked  
# with java -jar snpEff.jar databases
```

```
name: GRCh37.75
```

```
# Path to the reference genome, ideally as it is provided by the GATK bundle.
```

```
genome: /Users/patient/Desktop/hg19/hg19.fa
```

```
# Path to any database of known variants, ideally as it is provided by the GATK bundle.
```

```
known-variants: /Users/patient/Desktop/hg19/dbsnp_150.hg19.vcf.gz
```


Edit configuration files

config.yaml

filtering:

Set to true in order to apply machine learning based recalibration of
quality scores instead of hard filtering.

vqsr: false

hard:

hard filtering as outlined in GATK docs

(<https://gatkforums.broadinstitute.org/gatk/discussion/2806/howto-apply-hard-filters-to-a-call-set>)

snvs:

"QD < 2.0 || QUAL < 100.0 || DP < 50.0 || SOR > 3.0 || FS > 60.0 || MQ < 40.0 || MQRankSum < -12.5 ||
ReadPosRankSum < -8.0"

indels:

"QD < 2.0 || QUAL < 100 || DP < 50.0 || FS > 200.0 || ReadPosRankSum < -20.0"

#depth of coverage threshold to apply to variants identified with MuTect2

depth: "DP < 30"

Edit configuration files

config.yaml

```
processing:  
  remove-duplicates: true  
  # Uncomment and point to a bed file with, e.g., captured regions if necessary,  
  # see https://gatkforums.broadinstitute.org/gatk/discussion/4133/when-should-i-use-l-to-pass-in-a-list-of-intervals.  
  restrict-regions: /Users/patient/Desktop/bedfiles/CCP.20131001.designed.bed  
  # If regions are restricted, uncomment this to enlarge them by the given value in order to include  
  # flanking areas.  
  region-padding: 100
```

Edit configuration files

config.yaml

```
params:
```

```
...
```

```
trimmomatic:
```

```
  pe:
```

```
    trimmer:
```

```
    # See trimmomatic manual for adding additional options, e.g. for adapter trimming.
```

```
      - "LEADING:0"
```

```
      - "TRAILING:0"
```

```
  se:
```

```
    trimmer:
```

```
    # See trimmomatic manual for adding additional options, e.g. for adapter trimming.
```

```
      - "LEADING:0"
```

```
      - "TRAILING:0"
```

Edit configuration files

config.yaml

```
resources:
```

```
...
```

```
  trim_reads:
```

```
    mem: 32000
```

```
    threads: 16
```

```
    walltime: 480
```

Execution

1. Activate the conda environment

```
$ conda activate snakemake
```

2. Check dry run

```
$ snakemake --use-conda -n --until merge_calls filter_mutect_2
```

```
Job stats:
job          count  min threads  max threads
-----
call_variants      24          1          1
combine_calls      24          1          1
compose_regions    24          1          1
filter_mutect_2     1          1          1
filter_mutect_calls 1          1          1
genotype_variants   24          1          1
hard_filter_calls    2          1          1
map_reads           1          1          1
mark_duplicates      1          1          1
merge_calls         1          1          1
merge_variants      1          1          1
mutect              1          1          1
obtain_recal_table   1          1          1
recalibrate_base_qualities 1          1          1
samtools_index       1          1          1
select_calls        2          1          1
sort_bed            1          1          1
trim_reads_pe        1          1          1
total              112          1          1

This was a dry-run (flag -n). The order of jobs does not reflect the order of execution.
```

Execution

3. Launch Varca

```
$ snakemake --use-conda --jobs 4 --until merge_calls filter_mutect_2
```

It takes a while, so let's move to the other exercise while it finishes!